

## Retroviral generation of stable cell lines

### I. Generate retrovirus

Day 1

1. Transfect 293T cells (40-60% confluent) with 3 plasmids: 9 µg of retroviral plasmid (e.g. pBABE-ING2), 4.5 µg of GAG-POL plasmid and 4.5 µg of VSV plasmid using Transit293 transfection reagent.

Day 2

2. (Optional) The next day, change the medium, but only replace 4-5 ml DMEM-10% FCS (to keep virus concentrated). The medium change is required to remove a growth inhibitory factor produced by 293 cells upon transfection.

Day 3

3. 2 days after transfection, collect the supernatant (medium) with a 10 ml syringe, and add 4-5 ml fresh DMED.

Day 4

4. Filter the medium through **0.45µm** filter.
5. Aliquote 2 ml each and freeze @ -80oC or use fresh.
6. The next day, harvest the virus as day 3.

### II. Transduction

Day 3

1. Plate of HT1080 cells should be around 40-50% confluent (MMLV based retrovirus cannot enter non-dividing cells).
2. For a 10-cm dish, remove media and add 6 ml of new medium containing 64 µg Polybrene (16 µl of stock solution 4 mg/ml at -20oC).
3. Thaw the virus in the medium at 37oC quickly, and add 2ml into each plate (final polybrene conc. is 8 µg/ml).

Day 4

4. After 4-6 hrs, add 4 more mls media.
5. After 24 hrs, change media/and split if necessary.

### III. Stable line selection

Day 5

1. 24 hrs after the first media change (2 days after transduction), start selection with proper drugs (2 µg/ml puromycin (1 to 1000 dilution from 2 mg/ml stock) for HT1080, may vary for different cell types).

Day  
6-8

- (1) Split cells if necessary, cells need to be dividing for puro to be active.
- (2) Puro can be added 24 hrs or 48 hrs after transduction.
- (3) 1 mock transfected plate can be used to ensure that the puro did its job.
2. Change media every 1-2 days/split cells when necessary.
3. In 2-3 days all non-transduced cells should have been killed and you have your stable population of cells.